


Mini Centrifuge

Model: iSpin-7s



Multi-rotor, more working capacity



Brushless motor



High frequency and wide voltage



Low noise

TECHNICAL DATA

Model	iSpin-7s
Control mode	PI power direct drive
Max relative centrifugal force	3,286xg
Max rotation speed	7,000rpm±5%
Rotor	0.2/0.5/1.5/2ml*12 + 0.2ml*8*4 (Centrifuge tube) 0.5/1.0ml*12 (Micro blood tube)
Power	25W
Voltage	AC100-240V / DC24V

